LABQUALITY

External Quality Assessment Scheme

5-hydroxyindoleacetic Acid (5HIAA) Round 1, 2023

Specimens

Please find enclosed 2 liquid pooled human sera samples S001 and S002, 1 mL, each. No preservatives are added.

Caution

Quality control specimens derived from human blood must be handled with the same care as patient samples, i.e. as potential transmitters of serious diseases. The specimens are found to be HBsAg, HCVAb and HIVAgAb negative when tested with licensed reagents, but no known test method can offer complete assurance that the specimens will not transmit these or other infectious diseases.

Examinations

5-hydroxyindoleacetic Acid (5HIAA) (serum)

Storage and use

Samples are sent to participants as soluble, unfrozen ready-to-use preparations. Samples should be analyzed immediately upon arrival. If the assay cannot be performed immediately, samples should be frozen (-20 °C). Thaw samples shortly before analysis, allow to reach room temperature and mix until homogeneous. Analyze samples as patient samples.

Result reporting

Please enter the results and methods via LabScala (www.labscala.com). If you have any problems, please contact the EQA Coordinator.

Log in to LabScala, fill in the fields for sample arrival date and storage conditions. Then use the green button to enter the results page. Select your method from the drop-down menu. Please report your result as nmol/L. If your device gives the result as ng/mL, please convert the results to nmol/L using a conversion factor 5.231.

S001



S002



2023-04-18

INSTRUCTIONS

Product no. 2525 LQ751622011-012/FI Subcontracting: Sample preparation, sample pretesting

If the kit is incomplete or contains damaged specimens, please report immediately to info@labquality.fi.

The results should be reported no later than **May 15, 2023**.

Inquiries

EQA Coordinator Päivi Ranta paivi.ranta@labquality.fi

Labquality Oy

Kumpulantie 15 FI-00520 HELSINKI Finland

Tel. + 358 9 8566 8200 Fax + 358 9 8566 8280


info@labquality.fi www.labquality.com






LABQUALITY

Hydroxyindolyl acetate (5-), serum |LC-MS/MS






| | x _{pt} | sd | SEM | CV% | n |
|-------------|-----------------|-------|-------|------|---|
| LC-MS/MS | 270.01 nmol/l | 27.35 | 12.23 | 10.1 | 5 |
| All methods | 270.69 nmol/l | 24.52 | 10.01 | 9.1 | 6 |

| | x _{pt} | sd | SEM | CV% | n |
|-------------|-----------------|-------|-------|-----|---|
| LC-MS/MS | 560.10 nmol/l | 42.62 | 19.06 | 7.6 | 5 |
| All methods | 546.75 nmol/l | 50.23 | 20.51 | 9.2 | 6 |

History



▲ diff% ▽ Z-score is uncertain due to the small number of observations

| Round | Sample | ^X pt | Result | diff% | z-score |
|-------|---------------------|-----------------|---------|--------|---------|
| 23/1 | Sample S002-2525 | 560.10 | 611.00 | 9.09% | 1.19 |
| 23/1 | Sample S001-2525 | 270.01 | 288.00 | 6.66% | 0.66 |
| 22/2 | Sample S002-2525 | 100.75 | 176.60 | 75.29% | 2.01 |
| 22/2 | Sample S001-2525 | 3426.21 | 3395.50 | -0.90% | -0.11 |
| 22/1 | Sample S002-2525 | 539.56 | 580.90 | 7.66% | 0.90 |
| 22/1 | Sample S001-2525 | 94.30 | 108.50 | 15.05% | 1.03 |
| 21/2 | Sample S002-2525 | 149.51 | 174.54 | 16.74% | 0.85 |
| 21/2 | Sample S001-2525 | 169.67 | 187.59 | 10.56% | 0.49 |
| 21/1 | Sample S002-2525 | 232.09 | 264.80 | 14.09% | 0.74 |
| 21/1 | Sample S001-2525 | 568.71 | 569.10 | 0.07% | 0.00 |

17.05.2023



5-hydroxyindoleacetic Acid (5-HIAA), April, 1-2023



Report info

Participants

7 participants from 6 countries.

Report info

Your own result should be compared to others using the same method. Assigned values (x_{pt}, target values) are means of the results where results deviating more than +/- 3*standard deviation from the median are removed. The standard uncertainty (u) of

the assigned value is reported as standard error of the mean (SEM). Additionally, if the measurement uncertainty of the target value is large an automatic text is printed on the report: "The uncertainty of the assigned value is not negligible, and evaluations could be affected."
In case the client's result is the only one in the method group, no assigned value will be calculated, no target area shown, and no statistics calculated. In case there are only a few results in the client's own method group, the result can be compared to all method mean or to a group that is similar to the own method. Results reported with < or > -signs cannot be included in the statistics.


For information on report interpretation and performance evaluation, please see the "EQAS Interpretation guidelines" LabScala User instructions (top right corner? Help link).




Sample S001 | Hydroxyindolyl acetate (5-), serum, nmol/l

| Methodics | x _{pt} | Median | sd | CV% | SEM | min | max | Outliers | n |
|----------------------------------|-----------------|--------|-------|------|-------|--------|--------|----------|---|
| LC-MS/MS | 270.01 | 287.71 | 27.35 | 10.1 | 12.23 | 240.00 | 294.00 | - | 5 |
| Other, specify in comments field | - | - | - | - | - | 274.10 | 274.10 | - | 1 |
| All | 270.69 | 280.91 | 24.52 | 9.1 | 10.01 | 240.00 | 294.00 | _ | 6 |

Sample S001 | Hydroxyindolyl acetate (5-), serum, nmol/l| histogram summaries in LabScala









Sample S002 | Hydroxyindolyl acetate (5-), serum, nmol/l

| Methodics | x _{pt} | Median | sd | CV% | SEM | min | max | Outliers | n |
|----------------------------------|-----------------|--------|-------|-----|-------|--------|--------|----------|---|
| LC-MS/MS | 560.10 | 564.00 | 42.62 | 7.6 | 19.06 | 514.00 | 611.00 | - | 5 |
| Other, specify in comments field | - | _ | - | - | - | 480.00 | 480.00 | - | 1 |
| All | 546.75 | 542.21 | 50.23 | 9.2 | 20.51 | 480.00 | 611.00 | - | 6 |

Sample S002 | Hydroxyindolyl acetate (5-), serum, nmol/l| histogram summaries in LabScala





16.05.2023



5-hydroxyindoleacetic Acid (5-HIAA), April, 1-2023

Report info

Participants

7 participants from 6 countries.

Report info

Your own result should be compared to others using the same method. Assigned values (x_{pt}, target values) are means of the results where results deviating more than +/- 3*standard deviation from the median are removed. The standard uncertainty (u) of

the assigned value is reported as standard error of the mean (SEM). Additionally, if the measurement uncertainty of the target value is large an automatic text is printed on the report: "The uncertainty of the assigned value is not negligible, and evaluations could be affected."
In case the client's result is the only one in the method group, no assigned value will be calculated, no target area shown, and no statistics calculated. In case there are only a few results in the client's own method group, the result can be compared to all method mean or to a group that is similar to the own method. Results reported with < or > -signs cannot be included in the statistics.

For information on report interpretation and performance evaluation, please see the "EQAS Interpretation guidelines" LabScala User instructions (top right corner? Help link).

16.05.2023 3/3

LABQUALITY

External Quality Assessment Scheme

5-hydroxyindoleacetic Acid (5HIAA) Round 1, 2023

Specimens

Sample S001 (LQ751623011) and sample S002 (LQ751623012) were prepared by dilution of serum pools with elevated 5HIAA concentration into a serum matrix.

Based on the previous tests and the results of this round, the samples are homogeneous, stable, and suitable for the external quality assessment scheme.

The materials were sent without temperature control packaging.

Report info

Please see the description of the data analysis on the last page of the laboratory-specific histogram and Global report. It is important to read the Final report first, because it contains important information of the samples and results in each round.

Comments - EQA Coordinator

There were six participants on the round of which five used mass spectrometry and one HPLC gradient method. If there is only one participant in any method group the target area can't be given. The client was given a comment what is their difference % compared to the mean of all results.

Comments - Expert

Both serum samples of this round were prepared to contain an elevated concentration of 5HIAA. All laboratories performed well with results within the ± 15 % target range.

End of report

2023-06-22

FINAL REPORT

Product no. 2525

Subcontracting: Sample preparation, sample pretesting

 Samples sent
 2023-04-18

 Round closed
 2023-05-15

 Final report
 2023-06-22

Request for correction

Typing errors in laboratory's result forms are on laboratory's responsibility.

Labquality accepts responsibility only for result processing. Requests must be notified by writing within three weeks from the date of this letter.

Authorized by

EQA Coordinator Päivi Ranta paivi.ranta@labquality.fi

Expert

PhD Outi Itkonen HUS, Diagnostic Center, Helsinki, Finland.

Labquality Oy

Kumpulantie 15 FI-00520 HELSINKI Finland

Tel. + 358 9 8566 8200

info@labquality.fi www.labquality.com



Copyright © Labquality Oy

Labquality does not permit any reproduction for commercial purposes of any portion of the material subject to this copyright. Labquality prohibits any use of its name, or reference to Labquality EQA program, or material in this report in any advertising, brochures or other commercial publications. Labquality EQA data do not necessarily indicate the superiority of instruments, reagents, testing equipment or materials used by participating laboratories. Use of Labquality EQA data to suggest superiority or inferiority of equipment or materials may be deceptive and misleading. Proficiency test results are handled confidentially. Labquality will not issue any statements to third parties of the performance of laboratories in external quality assessment schemes unless otherwise agreed.