

# LABQUALITY

External Quality Assessment Scheme

## Blood culture, screening, round 1, 2018

Welcome to participate in the EQA scheme for Blood culture, screening.

### Specimens

Please find enclosed two lyophilized specimens and vials of rehydration fluid. Handle the specimens with the same care as corresponding clinical specimens capable of transmitting infectious disease. Specimens should after arrival be stored at 2 ... 8 °C. Follow the standard operating procedure of your laboratory for disposal of the specimens. Please follow the instructions, incubate and culture the specimens and read the results. Record your results and the methods used in the enclosed result forms.

### Background information

**Specimen 001** (S001: LQ761818011): Sepsis. Hospitalized after a car crash accident in Greece.

**Specimen 002** (S002: LQ761818012): Elderly patient with pneumonia.

### Handling instructions:

1. Let the specimen and the rehydration fluid warm up to room temperature.
2. Cut the foil packet open at the end where you can feel the thicker part of the loop.
3. Remove the plastic sheath from the loop. Break the loop shaft off from handle directly into the tube containing warm rehydration fluid (blue cap).
4. Incubate the tube for 30 minutes at 35 ... 37 °C.
5. Check that the black film dissolves completely out of the loop. Mix well and discard the loop.
6. Measure 10 mL of blood (taken from a healthy person or animal) into a sterile tube. \*
7. Add **10 µL** of bacterial specimen to the blood. \*\*
8. Mix well the content of the tube.
9. Divide the content of the tube into blood culture bottles: 5 mL in aerobic bottle and 5 mL in anaerobic bottle, or, if only one bottle is used (e.g. Oxoid Signal), add the whole content of the tube into this bottle.

### Please note:

\* To minimize the risk of coagulation, the blood can alternatively be added directly into the blood culture bottles as follows: 5 mL in aerobic bottle and 5 mL in anaerobic bottle, or, if only one bottle is used (e.g. Oxoid Signal) 10 mL directly into the bottle.

\*\* If the blood is added directly into the blood culture bottles proceed as follows; **add 10 µL of bacterial specimen** into 500 µL of 0.9% NaCl, mix well and add 250 µL of this bacterial specimen into each blood culture bottle, or if only one bottle is used, add all 500 µL into the bottle.

Incubate similar to clinical specimens.

2018-03-06

### INSTRUCTIONS

Product no. 5101

UN3373

LQ761818011-012/US

Subcontracting: sample pretesting

### The shipment includes

- 2 lyophilized specimens
- 2 vials of rehydration fluids
- 2 result forms

If the kit is incomplete or contains damaged specimens, please report immediately to the EQA coordinator in order to obtain replacements.

### Closing date

The results should be in Labquality not later than **April 3, 2018.**

### Expected results

The expected results of the round are displayed in LabSca-la in the "View reports" section on April 5, 2018.

### Inquires

EQA Coordinator

Yvonne Björkman

T. +358 50 501 4335

yvonne.bjorkman@labquality.fi

### Labquality

Kumpulantie 15

FI-00520 Helsinki FINLAND

Telephone

+358 9 8566 8200

Fax

+358 9 8566 8280

+358 9 8566 8281

info@labquality.fi

www.labquality.fi

© Labquality

**FINAS**  
Finnish Accreditation Service  
PT02 (EN ISO/IEC 17043)



**Reporting your results:**

The results should be recorded on the enclosed paper forms; electronic LabScala forms are not yet available. We recommend returning your result forms as an attachment to [info@labquality.fi](mailto:info@labquality.fi). Also ordinary mail or fax can be used. Detailed contact info in column on the reverse page.

**Filling in the paper form:**

Please mark the name and the client code of your laboratory in the top right hand corner of the result form.

Report the blood culture system used in your laboratory. Next report the length of incubation when growth was detected, or, bottles were discarded as negative. The result of Gram staining is recorded in the section "further handling". Finally tick the appropriate "report to the clinician". If your laboratory represents a stem cell bank and does not identify the growth detected, kindly tick the appropriate box.

**Closing date** of the round is shown in the column on the right side of reverse page. Kindly note, that results, which are received at Labquality after this date, are not accepted in the result processing.

**The expected results** of the round are displayed in LabScala in the "View reports" -section on the date shown in the column on the right side of reverse page.

All comments concerning the scheme in general or the specimens are most welcome.

**Barcodes for the specimens:**

S001: LQ761818011



S002: LQ761818012



## Blood culture, round 1, 2018

Welcome to participate in the EQA Scheme for Blood culture.

### Specimens

Please find enclosed two lyophilized specimens and vials of rehydration fluid. Handle the specimens with the same care as corresponding clinical specimens capable of transmitting infectious disease. Specimens should be stored after arrival at 2 ... 8 °C. Follow the standard operating procedure of your laboratory for disposal of the specimens. Please follow the instructions, incubate and culture the specimens and read the results. Record your results and the methods used in the enclosed result forms.

### Background information

**Specimen 001** (S001: LQ761818011): Sepsis. Hospitalized after a car crash accident in Greece.

**Specimen 002** (S002: LQ761818012): Elderly patient with pneumonia.

### Handling instructions:

1. Let the specimen and the rehydration fluid warm up to room temperature.
2. Cut the foil packet open at the end where you can feel the thicker part of the loop.
3. Remove the plastic sheath from the loop. Break the loop shaft off from handle directly into the tube containing warm rehydration fluid (blue cap).
4. Incubate the tube for 30 minutes at 35 ... 37 °C.
5. Check that the black film dissolves completely out of the loop. Mix well and discard the loop.
6. Measure 10 mL of blood (taken from a healthy person or animal) into a sterile tube.\*
7. Add **10 µL** of bacterial specimen to the blood.\*\*
8. Mix well the content of the tube.
9. Divide the content of the tube into blood culture bottles: 5 mL in aerobic bottle and 5 mL in anaerobic bottle, or, if only one bottle is used (e.g. Oxoid Signal), add the whole content of the tube into this bottle.

### Please note:

\* To minimize the risk of coagulation, the blood can alternatively be added directly into the blood culture bottles as follows: 5 mL in aerobic bottle and 5 mL in anaerobic bottle, or, if only one bottle is used (e.g. Oxoid Signal) 10 mL directly into the bottle.

\*\* If the blood is added directly into the blood culture bottles proceed as follows; **add 10 µL of bacterial specimen** into 500 µL of 0.9% NaCl, mix well and add 250 µL of this bacterial specimen into each blood culture bottle, or if only one bottle is used, add all 500 µL into the bottle.

Incubate similar to clinical specimens.

2018-03-06

### INSTRUCTIONS

Product no. 5100  
UN3373  
LQ761818011-012/US  
Subcontracting: sample pretesting

### The shipment includes

- 2 lyophilized specimens
- 2 vials of rehydration fluids.
- result forms and a code list

If the kit is incomplete or contains damaged specimens, please report immediately to the EQA coordinator in order to obtain replacements.

### Closing date

The results should be in Labquality not later than **April 3, 2018.**

### Expected results

The expected results of the round are displayed in LabSca-la in the "View reports" section on April 5, 2018.

### Inquiries

EQA Coordinator  
Yvonne Björkman  
yvonne.bjorkman@labquality.fi  
T. +358 50 501 4335

### Labquality

Kumpulantie 15  
FI-00520 Helsinki FINLAND

### Telephone

+358 9 8566 8200

### Fax

+358 9 8566 8280

+358 9 8566 8281

info@labquality.fi  
www.labquality.fi

© Labquality

  
**FINAS**  
Finnish Accreditation Service  
PT02 (EN ISO/IEC 17043)



**Reporting your results:**

The results should be recorded on the enclosed paper forms; electronic LabScala forms are not yet available. We recommend returning your result forms as an attachment to [info@labquality.fi](mailto:info@labquality.fi). Also ordinary mail or fax can be used. Detailed contact info in column on the reverse page.

**Filling in the paper form:**

Please mark the name and the client code of your laboratory in the top right hand corner of the result form.

Report the blood culture system used in your laboratory. Next report the length of incubation when growth was detected, or, bottles were discarded as negative.

In the section "further handling", please report the results of gram staining.

The final report should include the answer in written and the corresponding code number found on the enclosed code list. If the answer is not found on the list, please use code 9999.

**Susceptibility testing (only for specimen 001):**

In order to be comparable with the other participants' results, please record which standard is followed in your laboratory for susceptibility testing procedures.

Susceptibility test results are given only for the antimicrobial agents that are routinely used in your laboratory for the isolated microbe in question.

For the disk diffusion method, report the inhibitory zone diameter (mm) and if MIC method is used, report the MIC value as mg/L. In the last column report the corresponding SIR-interpretation (S/I/R). The interpretation should be reported by taking into consideration the possible resistance mechanisms of the microbe.

**Closing date** of the round is shown in the column on the right side of reverse page. Kindly note, that results, which are received at Labquality after this date, are not accepted in the result processing.

**The expected results** of the round are published in LabScala in the "View reports" -section on the date shown in the column on the right side of reverse page.

**Specimen barcodes:**

S001: LQ761818011



S002: LQ761818012



When you have identified the microbe or obtained the gram stain result, please select the corresponding code number from this list. Transfer the number to the result form. If your identification of the microbe or the gram staining result can not be found on this code list, please use code 9999 and write the name of the microbe / staining result.

CODE NAME	CODE NAME	CODE NAME
8556 <i>Achromobacter xylosoxidans</i> ssp. <i>denitrificans</i>	7640 <i>Escherichia coli</i>	8490 <i>Shigella doydii</i>
8038 <i>Acinetobacter</i> sp.	5688 <i>Escherichia coli</i> , EHEC	3667 <i>Shigella sonnei</i>
7647 <i>Acinetobacter baumannii</i>	7479 <i>Finegoldia magna</i> (syn. <i>Peptostreptococcus magnus</i> )	7428 <i>Sphingobacterium multivorum</i>
5873 <i>Acinetobacter calcoaceticus</i>	1309 <i>Flavobacterium</i> sp.	9086 <i>Staphylococcus</i> sp.
2428 <i>Acinetobacter lwoffii</i>	2092 <i>Fusobacterium</i> sp.	2698 <i>Staphylococcus aureus</i>
3526 <i>Actinobacillus actinomycetemcomitans</i>	3393 <i>Fusobacterium necrophorum</i>	4664 <i>Staphylococcus capitis</i>
4845 <i>Actinomyces israelii</i>	9080 <i>Fusobacterium nucleatum</i>	9015 <i>Staphylococcus epidermidis</i>
3421 <i>Actinomyces odontolyticus</i>	8491 <i>Gardnerella vaginalis</i>	6193 <i>Staphylococcus haemolyticus</i>
8795 <i>Actinomyces</i> sp.	8409 <i>Gemella morbillorum</i>	5736 <i>Staphylococcus hominis</i>
3396 <i>Aerococcus</i> sp.	0998 <i>Haemophilus</i> sp.	7270 <i>Staphylococcus lugdunensis</i>
3354 <i>Aeromonas caviae</i>	7572 <i>Haemophilus influenzae</i>	3889 <i>Staphylococcus saprophyticus</i>
5837 <i>Aeromonas</i> sp.	5161 <i>Haemophilus influenzae</i> , type b	8249 <i>Staphylococcus warneri</i>
4301 <i>Aeromonas hydrophila</i>	5505 <i>Haemophilus parainfluenzae</i>	8805 <i>Staphylococcus xylosus</i>
9911 <i>Alcaligenes</i> sp.	5316 <i>Kingella kingae</i>	1885 <i>Stenotrophomonas maltophilia</i>
8123 <i>Alcaligenes faecalis</i>	9597 <i>Klebsiella</i> sp.	1933 <i>Streptococcus</i> sp.
7694 <i>Arcanobacterium haemolyticum</i>	9416 <i>Klebsiella oxytoca</i>	1933-4 <i>Streptococcus</i> sp. (β-hem., not Group A)
7403 <i>Aspergillus</i> sp.	1916 <i>Klebsiella ozaenae</i>	1769 <i>Streptococcus agalactiae</i> (β-hem., Group B)
1825 <i>Aspergillus fumigatus</i>	8712 <i>Klebsiella pneumoniae</i>	5158-1 <i>Streptococcus anginosus</i> -group
1415 <i>Bacillus</i> sp.	6393 <i>Lactobacillus acidophilus</i>	7337 <i>Streptococcus anginosus</i>
2216 <i>Bacillus cereus</i>	3498 <i>Leuconostoc</i> sp.	9074 <i>Streptococcus bovis</i>
8738 <i>Bacillus subtilis</i>	2963 <i>Listeria</i> sp.	1018 <i>Streptococcus</i> sp., β-hem., Group C
5938 <i>Bacteroides</i> sp.	9370 <i>Listeria monocytogenes</i>	1053 <i>Streptococcus</i> sp., β-hem., Group G
8059 <i>Bacteroides fragilis</i>	6991 <i>Micrococcus</i> sp.	5158-1 <i>Streptococcus milleri</i> -group
8065 <i>Bacteroides ovatus</i>	8906 <i>Moraxella</i> sp.	9327 <i>Streptococcus mitis</i> -group
1136 <i>Bacteroides uniformis</i>	9463 <i>Moraxella catarrhalis</i>	7986 <i>Streptococcus mitis</i>
6954 <i>Bacteroides vulgatus</i>	8921 <i>Moraxella osloensis</i>	3278 <i>Streptococcus mutans</i>
9972 <i>Bifidobacterium</i> sp.	2249 <i>Morganella morganii</i>	6147 <i>Streptococcus pneumoniae</i>
7071 <i>Bordetella parapertussis</i>	6485 <i>Mycobacterium</i> sp.	4030 <i>Streptococcus pyogenes</i> (β-hem. Group A)
5942 <i>Bordetella pertussis</i>	1026 <i>Mycobacterium abscessus</i>	7477 <i>Streptococcus salivarius</i>
6415 <i>Brevundimonas diminuta</i>	3342 <i>Mycobacterium tuberculosis</i>	7738 <i>Streptococcus sanguinis</i>
1637 <i>Burkholderia cepacia</i>	9640 <i>Neisseria</i> sp.	9504 <i>Streptococcus viridans</i> -group
8847 <i>Campylobacter</i> sp.	1757 <i>Neisseria gonorrhoeae</i>	5953 <i>Veillonella</i> sp.
6658 <i>Campylobacter coli</i>	5335 <i>Neisseria lactamica</i>	4751 <i>Veillonella parvula</i>
5342 <i>Campylobacter jejuni</i>	7202 <i>Neisseria meningitidis</i>	9434 <i>Vibrio</i> sp.
1119 <i>Candida albicans</i>	9994 <i>Neisseria mucosa</i>	1755 <i>Vibrio cholerae</i>
3034 <i>Candida glabrata</i>	5022 <i>Nocardia</i> sp.	6827 Yeast, other than <i>Candida albicans</i>
2828 <i>Candida krusei</i>	7060 <i>Nocardia asteroides</i>	3459 <i>Yersinia</i> sp.
6847 <i>Candida</i> sp.	1156 <i>Ochrobactrum anthropi</i>	8402 <i>Yersinia enterocolitica</i>
1013 <i>Capnocytophaga canimorsus</i>	8310 <i>Oligella ureolytica</i>	8258 <i>Yersinia pseudotuberculosis</i>
7665 <i>Cardiobacterium hominis</i>	7260 <i>Oligella urethralis</i>	
8975 <i>Citrobacter</i> sp.	5214 <i>Parvimonas micra</i>	
6168 <i>Citrobacter diversus</i>	7008 <i>Pasteurella</i> sp.	
3517 <i>Citrobacter freundii</i>	6536 <i>Pasteurella multocida</i>	
4933 <i>Clostridium</i> sp.	5507 <i>Peptostreptococcus</i> sp.	
9607 <i>Clostridium difficile</i>	5494 <i>Peptostreptococcus anaerobius</i>	
1421 <i>Clostridium innocuum</i>	6261 <i>Plesiomonas shigelloides</i>	
9793 <i>Clostridium perfringens</i>	8227 <i>Prevotella melaninogenica</i>	
3359 <i>Clostridium</i> sp.	2913 <i>Propionibacterium</i> sp.	
6715 <i>Clostridium septicum</i>	6289 <i>Propionibacterium acnes</i>	
7695 <i>Clostridium sporogenes</i>	2347 <i>Proteus</i> sp.	
4797 <i>Clostridium tertium</i>	6986 <i>Proteus mirabilis</i>	
8334 <i>Corynebacterium</i> sp.	9831 <i>Proteus vulgaris</i>	
4405 <i>Corynebacterium diphtheriae</i>	6283 <i>Providencia</i> sp.	
1190 <i>Corynebacterium jeikeium</i>	6872 <i>Providencia stuartii</i>	
6499 <i>Corynebacterium urealyticum</i>	8289 <i>Pseudomonas</i> sp.	
6064 <i>Cronobacter sakazakii</i> (syn. <i>Enterobacter sakazakii</i> )	7093 <i>Pseudomonas aeruginosa</i>	
8068 <i>Cryptococcus neoformans</i>	9574 <i>Pseudomonas fluorescens</i>	
6538 <i>Eggerthella lenta</i> (syn. <i>Eubacterium lentum</i> )	5462 <i>Pseudomonas putida</i>	
9293 <i>Eikenella corrodens</i>	9340 <i>Ralstonia pickettii</i> (syn. <i>Pseudomonas pickettii</i> )	
1605 <i>Enterobacter</i> sp.	7800 <i>Rhodococcus equi</i>	
4596 <i>Enterobacter aerogenes</i>	8672 <i>Rothia mucilaginosa</i> (syn. <i>Stomatococcus mucilaginosus</i> )	
4937 <i>Enterobacter agglomerans</i> -group	8019 <i>Salmonella</i> sp.	
1841 <i>Enterobacter cloacae</i>	6015 <i>Salmonella</i> Enteritidis	
5696 <i>Enterococcus</i> sp.	5887 <i>Salmonella</i> Typhimurium	
4547 <i>Enterococcus casseliflavus</i>	8244 <i>Serratia</i> sp.	
8864 <i>Enterococcus faecalis</i>	7155 <i>Serratia marcescens</i>	
3093 <i>Enterococcus faecium</i>	7058 <i>Serratia rubidaea</i>	
3129 <i>Enterococcus gallinarum</i>	1860 <i>Shewanella putrefaciens</i>	
4712 <i>Erysipelothrix rhusiopathiae</i>	4468 <i>Shigella</i> sp.	
	2832 <i>Shigella dysenteriae</i>	
	5614 <i>Shigella flexneri</i>	
		<b>PRELIMINARY IDENTIFICATION: (G:327)</b>
		G010 Aerobe grampositive cocci
		G011 Aerobe grampositive cocci in chains
		G012 Aerobe grampositive diplococci
		G013 Aerobe grampositive cocci in clusters
		G016 Aerobe grampositive bacilli
		G017 Aerobe grampositive sporeforming bacilli
		G020 Aerobe gramnegative bacilli
		G021 Aerobe gramnegative diplococci
		G022 Aerobe gramnegative vibrio
		G025 Anaerobe grampositive cocci
		G026 Anaerobe grampositive bacilli
		G027 Anaerobe grampositive sporeforming bacilli
		G028 Anaerobe gramnegative cocci
		G029 Anaerobe gramnegative bacilli
		0320 Gramnegative bacteria
		0321 Grampositive bacteria
		6369 Yeast
		0341 Mould
		<b>OTHER:</b>
		0331 Identification not performed in this laboratory, referred to another laboratory
		0333 Normal flora
		0334 Mixed flora
		0337 No growth (sterile)
		0338 No final result
		9999 Other; please specify on the result form