

#### Specimen S001 | Calprotectin, quantitative, mg/kg

Methodics	<sup>X</sup> pt	Median	sd	CV%	SEM	min	max	Outliers	n
Accysay Diagnostics, CalproGo (50 mg/kg)	-	-	-	-	-	6.0	6.0	-	1
BIOHIT Calprotectin Elisa (50 mg/kg)	-	-	-	-	-	28.3	28.3	-	1
Bühlmann fCal ELISA (50 mg/kg)	31.9	31.9	11.2	35.0	7.9	24.0	39.8	-	2
Bühlmann fCal turbo (50 mg/kg)	30.0	25.0	11.7	39.0	4.4	20.0	50.0	-	7
Bühlmann Quantum Blue fCal (50 mg/kg)	35.7	36.0	5.5	15.4	3.2	30.0	41.0	-	3
Calpro, CalproLab Elisa (ALP) (50 mg/kg)	29.7	29.7	6.4	21.4	4.5	25.2	34.2	-	2
Certest Biotec Turbilatex (50 mg/kg)	-	-	-	-	-	30.0	30.0	-	1
DiaSorin, Liaison Clia (50 mg/kg)	18.4	18.4	4.9	26.5	1.4	11.0	27.1	-	13
Euroimmun, ELISA (50 mg/kg)	37.0	37.0	2.8	7.6	2.0	35.0	39.0	-	2
Immundiagnostik, PhiCal ELISA (50 mg/kg)	51.5	51.5	0.7	1.4	0.5	51.0	52.0	-	2
Inova Diagnostics, QUANTA Lite (50 mg/kg)	19.8	19.8	1.4	7.1	1.0	18.8	20.8	-	2
Orgentec, ELISA (50 mg/kg)	18.5	18.5	5.1	27.5	3.6	14.9	22.1	-	2
Sentinel Diagnostics, CALiaGold (50 mg/kg)	-	-	-	-	-	26.6	26.6	-	1
Thermo Fisher, Phadia Elia 2 (50 mg/kg)	16.8	16.0	4.9	29.1	1.5	11.0	27.0	-	10
All	24.3	21.6	10.9	44.8	1.6	6.0	52.0	-	49

18.06.2018



#### Specimen S001 | Calprotectin, qualitative, -

Methodics	Negative	Total
Certest Biotec (50 mg/kg)	2	2
Nal von Minden, Nadal calprotectin (50 mg/kg)	2	2
VIDIA spol.r.s.o Rapid-VIDITEST (50 mg/kg)	1	. 1
Total	5	5

Copyright © Labquality Oy
Labquality does not permit any reproduction for commercial purposes of any portion of the material subject to this copyright. Labquality prohibits any use of its name, or reference to Labquality EQA program, or material in this report in any advertising, brochures or other commercial publications. Labquality EQA data do not necessarily indicate the superiority of instruments, reagents, testing equipments or materials used by participating laboratories. Use of Labquality EQA data to suggest superiority or inferiority of equipments or materials may be deceptive and misleading.

Proficiency test results are handled confidentially. Labquality will not issue any statements to third parties of the performance of laboratories in external quality assessment schemes unless otherwise agreed.

2/4 18.06.2018



#### Specimen S002 | Calprotectin, quantitative, mg/kg

Methodics	<sup>x</sup> pt	Median	sd	CV%	SEM	min	max	Outliers	n
Accysay Diagnostics, CalproGo (50 mg/kg)	-	-	-	-	-	42.0	42.0	-	1
BIOHIT Calprotectin Elisa (50 mg/kg)	-	-	-	-	-	147.4	147.4	-	1
Bühlmann fCal ELISA (50 mg/kg)	204.6	204.6	87.1	42.6	61.6	143.0	266.2	-	2
Bühlmann fCal turbo (50 mg/kg)	146.3	148.8	29.0	19.8	11.8	110.0	191.0	-	6
Bühlmann Quantum Blue fCal (50 mg/kg)	117.3	96.5	57.9	49.4	29.0	74.0	202.0	-	4
Calpro, CalproLab Elisa (ALP) (50 mg/kg)	129.0	129.0	24.0	18.6	17.0	112.0	146.0	-	2
Certest Biotec Turbilatex (50 mg/kg)	-	-	-	-	-	267.0	267.0	-	1
DiaSorin, Liaison Clia (50 mg/kg)	111.2	102.0	22.0	19.8	6.1	87.0	155.0	-	13
Euroimmun, ELISA (50 mg/kg)	297.5	297.5	10.6	3.6	7.5	290.0	305.0	-	2
Eurospital, CalFast, ELISA (70 mg/kg)	-	-	-	-	-	105.0	105.0	-	1
Eurospital, Calprest, ELISA (50 mg/kg)	-	-	-	-	-	61.0	61.0	-	1
Immundiagnostik, PhiCal ELISA (50 mg/kg)	161.5	161.5	3.5	2.2	2.5	159.0	164.0	-	2
Inova Diagnostics, QUANTA Lite (50 mg/kg)	120.2	120.2	3.4	2.8	2.4	117.8	122.6	-	2
Orgentec, ELISA (50 mg/kg)	180.9	180.9	83.5	46.2	59.1	121.8	239.9	-	2
Sentinel Diagnostics, CALiaGold (50 mg/kg)	-	-	-	-	-	153.1	153.1	-	1
Thermo Fisher, Phadia Elia 2 (50 mg/kg)	151.7	157.0	41.0	27.0	13.0	84.0	210.0	-	10
All	139.3	132.5	52.4	37.6	7.4	42.0	290.0	1	51

18.06.2018



#### Specimen S002 | Calprotectin, qualitative, -

Methodics	Positive	Total
Certest Biotec (50 mg/kg)		2
Nal von Minden, Nadal calprotectin (50 mg/kg)		2
VIDIA spol.r.s.o Rapid-VIDITEST (50 mg/kg)		1
Total		5

Copyright © Labquality Oy
Labquality does not permit any reproduction for commercial purposes of any portion of the material subject to this copyright. Labquality prohibits any use of its name, or reference to Labquality EQA program, or material in this report in any advertising, brochures or other commercial publications. Labquality EQA data do not necessarily indicate the superiority of instruments, reagents, testing equipments or materials used by participating laboratories. Use of Labquality EQA data to suggest superiority or inferiority of equipments or materials may be deceptive and misleading.

Proficiency test results are handled confidentially. Labquality will not issue any statements to third parties of the performance of laboratories in external quality assessment schemes unless otherwise agreed.

18.06.2018 4/4