

## Antistreptolysin 1, 2018

### Specimens

The specimens **S001: LQ757118011** and **S002: LQ757118012** are to be analyzed for antistreptolysin and they are from adults. The specimens are originated from human and their volume is ca. 0.4 mL.

S001: LQ757118011



S002: LQ757118012



### Note

The specimens are HBsAg-, HCVAb and HIVAb negative according to the supplier. However, no known test method offers complete assurance that hepatitis B virus, HIV or other infectious agents are absent. Accordingly, these human blood-based products should be handled as potentially infectious human sera.

### Storage and use

The specimens should be stored at 2 ... 8 °C and analysed as soon as possible, preferably within a week.

**Please tap the tubes before opening to prevent loss of samples.**

### Results

Please report **your method information and results only** via the **electronic LabScala result form**, which is available at [www.labquality.fi](http://www.labquality.fi) (-> choose LabScala and login). The guidance for reporting the results electronically you will find also in LabScala behind the "LabScala User instructions" button.

In case you have problems entering your results, please immediately contact Labquality.

LabScala is easy to use:

1. Log in with your username and password.
2. Add your devices to My devices (if you haven't done that already).
3. Choose round and add your contact person(s).
4. Enter your results.  
If you want to reply more than one result, press the + (plus) button on the right hand. A new result field will appear.
5. Check your results.
6. Send results to Labquality.

2018-02-12

### INSTRUCTIONS

Product no. 5840  
LQ757118011-012/US, CH

**The shipment includes**  
-2 specimens

If the kit is incomplete or contains damaged specimens, please report immediately to the EQA Coordinator.

**Please note:**  
**Return your results via**  
**LabScala no later than**  
**March 9, 2018.**

The expected results of the round are published in LabScala in the View Reports section on March 14, 2018.

### Inquiries

EQA Coordinator  
Satu Eklund  
T. +358 45 773 107 87  
[satu.eklund@labquality.fi](mailto:satu.eklund@labquality.fi)

**Labquality**  
**Kumpulantie 15**  
**FI-00520 HELSINKI**

**Telephone**  
**+358 9 8566 8200**

**Fax**  
**+358 9 8566 8280**  
**+358 9 8566 8281**

[info@labquality.fi](mailto:info@labquality.fi)  
[www.labquality.fi](http://www.labquality.fi)

© Labquality

  
**FINAS**  
Finnish Accreditation Service  
PT02 (EN ISO/IEC 17043)

