


GLOBAL REPORT


	No of participants	No of responded participants	Response percentage
Helicobacter pylori, antibodies, December, 4-2017	59	52	88.1 %

Summary




Summary	AVR success rate
Sample S001	98.3 %
Sample S002	96.6 %
Average:	97.4 %

Sample S001




Sample S001 results	Responded	AVR success rate	Count
	Clinical interpretation, qualitative methods	100 %	12
	Clinical interpretation, quantitative methods	100 %	42
	HepyAb	100 %	12
	HepyAbA	70 %	10
	HepyAbG	100 %	42
Total:		98.3 %	118

Sample S001 Clinical interpretation, qualitative methods



Active or past infection

Sample S001 Clinical interpretation, quantitative methods




Active or past infection


Clinical interpretation, qualitative methods	Interpretation	Interpretation count	AVR success rate	Interpretation score
	Active or past infection	12	100 %	4
	Total:	12	100 %	

Clinical interpretation, quantitative methods	Interpretation	Interpretation count	AVR success rate	Interpretation score
	Active or past infection	42	100 %	4
	Total:	42	100 %	


Sample S001 HepyAb



Sample S001 HepyAbA



Sample S001 HepyAbG




HepyAb	Interpretation	Test	Interpretation count	Test count	AVR success rate	Interpretation score
	Positive		12		100 %	2
		AccuBiotech ACCU-TELL Rapid H. pylori Antibody Test		1		
		Biomaxima HELICO H. pylori Ab		1		
		Liaison H. pylori Ab (Diasorin)		2		
		Meridian Bioscience Immunocard H. pylori Ab		1		
		Meridian Bioscience Premier H. pylori Ab		5		
		Standard Diagnostics SD Bioline H. pylori Ab		1		
		SureScreen Diagnostics H. pylori Ab Cassette		1		
	Total:		12		100 %	

HepyAbA	Interpretation	Test	Interpretation count	Test count	AVR success rate	Interpretation score
	Negative		7		100 %	2
		Euroimmun Anti-Helicobacter pylori CagA IgA ELISA		1		
		Euroimmun Helicobacter pylori IgA EIA/ELISA		3		
		IBL Helicobacter pylori IgA ELISA		1		
		Siemens Enzygnost Anti-Helicobacter pylori IgA EIA		2		
	Borderline		1		0 %	0
		Euroimmun Helicobacter pylori IgA EIA/ELISA		1		
	Positive		2		0 %	0
		DiaPro Diagnostics HP IgA ELISA		1		
		Orgentec Alegria Anti-Helicobacter pylori IgA ELISA		1		
	Total:		10		70 %	


HepyAbG	Interpretation	Test	Interpretation count	Test count	AVR success rate	Interpretation score
	Positive		42		100 %	2
		Biohit GastroPanel H. pylori IgG ELISA		1		
		bioMerieux Vidas H. pylori IgG		5		
		DiaPro Diagnostics HP IgG ELISA		1		
		Diasorin Liaison H. pylori IgG		4		
		Diesse Chorus H. pylori IgG ELISA		1		
		Euroimmun Helicobacter pylori IgG EIA/ELISA		12		
		In-house H. pylori IgG ELISA		1		
		Inova Diagnostics QUANTA Lite Helicobacter pylori ELISA		2		
		Meridian Bioscience Premier H. pylori IgG EIA		1		
		Orgentec Alegria Anti-Helicobacter pylori IgG ELISA		1		
		SERION classic Helicobacter pylori IgG ELISA		1		
		Siemens Enzygnost Anti-Helicobacter pylori IgG EIA		6		
		Siemens Immulite H. pylori IgG		4		
		Vircell Helicobacter pylori VirClaia IgG		1		
		ZEUS H. pylori IgG ELISA		1		
	Total:		42		100 %	

Sample S002




Sample S002 results	Responded	AVR success rate	Count
	Clinical interpretation, qualitative methods	100 %	12
	Clinical interpretation, quantitative methods	96.4 %	42
	HepyAb	100 %	12
	HepyAbA	-	10
	HepyAbG	94 %	42
	Total:	96.6 %	118

Sample S002 Clinical interpretation, qualitative methods



Active or past infection

Sample S002 Clinical interpretation, quantitative methods




No indication of infection
Active or past infection
Borderline interpretation


Clinical interpretation, qualitative methods	Interpretation	Interpretation count	AVR success rate	Interpretation score
	Active or past infection	12	100 %	4
	Total:	12	100 %	

Clinical interpretation, quantitative methods	Interpretation	Interpretation count	AVR success rate	Interpretation score
	No indication of infection	1	0 %	0
	Borderline interpretation	1	50 %	2
	Active or past infection	40	100 %	4
	Total:	42	96.4 %	


Sample S002 HepyAb



Sample S002 HepyAbA



Sample S002 HepyAbG



HepyAb	Interpretation	Test	Interpretation count	Test count	AVR success rate	Interpretation score
	Positive		12		100 %	2
		AccuBiotech ACCU-TELL Rapid H. pylori Antibody Test		1		
		Biomaxima HELICO H. pylori Ab		1		
		Liaison H. pylori Ab (Diasorin)		2		
		Meridian Bioscience Immunocard H. pylori Ab		1		
		Meridian Bioscience Premier H. pylori Ab		5		
		Standard Diagnostics SD Bioline H. pylori Ab		1		
		SureScreen Diagnostics H. pylori Ab Cassette		1		
	Total:		12		100 %	

HepyAbA	Interpretation	Test	Interpretation count	Test count	AVR success rate	Interpretation score
	Negative		4		-	-
		Euroimmun Helicobacter pylori IgA EIA/ELISA		1		
		Orgentec Alegria Anti-Helicobacter pylori IgA ELISA		1		
		Siemens Enzygnost Anti-Helicobacter pylori IgA EIA		2		
	Positive		6		-	-
		DiaPro Diagnostics HP IgA ELISA		1		
		Euroimmun Anti-Helicobacter pylori CagA IgA ELISA		1		
		Euroimmun Helicobacter pylori IgA EIA/ELISA		3		
		IBL Helicobacter pylori IgA ELISA		1		
	Total:		10			

HepyAbG	Interpretation	Test	Interpretation count	Test count	AVR success rate	Interpretation score
	Negative		2		0 %	0
		DiaPro Diagnostics HP IgG ELISA		1		
		Euroimmun Helicobacter pylori IgG EIA/ELISA		1		
	Borderline		1		50 %	1
		Siemens Enzygnost Anti-Helicobacter pylori IgG EIA		1		
	Positive		39		100 %	2
		Biohit GastroPanel H. pylori IgG ELISA		1		
		bioMerieux Vidas H. pylori IgG		5		
		Diasorin Liaison H. pylori IgG		4		
		Diesse Chorus H. pylori IgG ELISA		1		
		Euroimmun Helicobacter pylori IgG EIA/ELISA		11		
		In-house H. pylori IgG ELISA		1		
		Inova Diagnostics QUANTA Lite Helicobacter pylori ELISA		2		
		Meridian Bioscience Premier H. pylori IgG EIA		1		
		Orgentec Alegria Anti-Helicobacter pylori IgG ELISA		1		
		SERION classic Helicobacter pylori IgG ELISA		1		
		Siemens Enzygnost Anti-Helicobacter pylori IgG EIA		5		
		Siemens Immulite H. pylori IgG		4		

Copyright © Labquality Oy
Labquality does not permit any reproduction for commercial purposes of any portion of the material subject to this copyright. Labquality prohibits any use of its name, or reference to Labquality EQA program, or material in this report in any advertising, brochures or other commercial publications. Labquality EQA data do not necessarily indicate the superiority of instruments, reagents, testing equipments or materials used by participating laboratories. Use of Labquality EQA data to suggest superiority or inferiority of equipments or materials may be deceptive and misleading. Proficiency test results are handled confidentially. Labquality will not issue any statements to third parties of the performance of laboratories in external quality assessment schemes unless otherwise agreed.

		Vircell Helicobacter pylori VirClia IgG		1		
		ZEUS H. pylori IgG ELISA		1		
	Total:		42		94 %	