

## INSTRUCTIONS

### *Borrelia burgdorferi* antibodies 1, 2018

#### Specimens

Enclosed you find two unknown human specimens  
S001: LQ757318011 and S002: LQ757318012, which are obtained from adults. The sample volume is ca. 0.5 mL.

S001: LQ757318011



S002: LQ757318012



#### Note

The specimens are HBs-Ag- and HIV-AgAb-negative according to the supplier. However, no known test method offers complete assurance that hepatitis B virus, HIV or other infectious agents are absent. Accordingly, these human blood-based products should be handled as potentially infectious human sera.

#### Storage and use

The specimens should be stored at 2 ... 8 °C and analysed as soon as possible, preferably within a week.

#### Results

Please report your **method information** and **results** only via LabScala.

The results should be recorded the same way as adult patient specimens. If you have any problems in reporting results via LabScala, please contact EQA Coordinator.

Product no. 5960  
LQ757318011-012/UK

This shipment includes  
2 specimens for analysis  
and instructions.

If the kit is incomplete or  
contains damaged speci-  
mens, please report im-  
mediately to the EQA Co-  
ordinator.

#### Please note:

The results should be  
reported via LabScala  
no later than  
**February 23, 2018.**

The expected results of the  
round are published in  
LabScala in the  
View Reports section on  
February 28, 2018.

#### Inquiries

EQA Coordinator  
Satu Eklund  
Tel. +358 45 773 107 87  
satu.eklund@labquality.fi

**Labquality**  
Kumpulantie 15  
FI-00520 HELSINKI

**Telephone**  
+358 9 8566 8200

**Fax**  
+358 9 8566 8280  
+358 9 8566 8281

info@labquality.fi  
www.labquality.fi

© Labquality

  
Finnish Accreditation Service  
PT02 (EN ISO/IEC 17043)

