

## *Bordetella pertussis* antibodies 1, 2018

### Specimens

Enclosed you find two unknown human serum specimens  
S001: LQ757218011 and S002: LQ757218012, which are obtained from  
adults. The specimen volume is ca. 0.3 mL.

### Note

The specimens are HBs-Ag- and HIV-AgAb-negative according to the  
supplier. However, no known test method offers complete assurance that  
hepatitis B virus, HIV or other infectious agents are absent. Accordingly,  
these human blood-based products should be handled as potentially in-  
fectious human sera.

### Storage and use

The specimens should be stored at 2 ... 8 °C and analysed as soon as  
possible, preferably within a week.

### Analytes

S- BopeAbA	S- Bordetella pertussis, IgA antibodies
S- BopeAbG	S- Bordetella pertussis, IgG antibodies
S- BopeAbM	S- Bordetella pertussis, IgM antibodies

### Results

Please report your **method information** and **results** only  
via **LabScala**.

The results should be recorded the same way as adult patient  
specimens.

### Barcodes



2018-01-29

### INSTRUCTIONS

Product no. 5950  
LQ757218011-012/UK

The round kit contains  
2 specimens for analysis  
and instructions.

If the kit is incomplete or  
contains damaged speci-  
mens, please report imme-  
diately to the EQA Coordi-  
nator.

### Please note:

The results should be reported  
via Labscala no later than  
**February 23, 2018.**

The expected results of the  
round are published in  
LabScala in the View Reports  
section on February 28, 2018.

### Inquiries

EQA Coordinator  
Satu Eklund  
Tel. +358 45 773 107 87  
satu.eklund@labquality.fi

**Labquality**  
Kumpulantie 15  
FI-00520 HELSINKI

**Telephone**  
+358 9 8566 8200

**Fax**  
+358 9 8566 8280  
+358 9 8566 8281

info@labquality.fi  
www.labquality.fi

© Labquality

